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Subject: Replace the products of Batch SE2010037, SE2011037

Date: Jan 06, 2020

### **The affected devices in global market:**

Name of device	Batch	Quantity	Catalog No.
VivaDiag SARS-CoV-2 Ag Rapid Tests	SE2010037	156,525 pcs	VCD05-01-011
VivaDiag SARS-CoV-2 Ag Rapid Tests	SE2011044	269,350 pcs	VCD05-01-011
VivaDiag SARS-CoV-2 Ag Rapid Tests	SE2011037	362,850 pcs	VCD05-01-011

### **The affected devices in Italy market:**

Name of device	Batch	Quantity	Catalog No.
VivaDiag SARS-CoV-2 Ag Rapid Tests	SE2010037	16,050 pcs	VCD05-01-011
VivaDiag SARS-CoV-2 Ag Rapid Tests	SE2011037	30,250 pcs	VCD05-01-011

### **Description of the problem:**

ANSM (case number: I2014389) has informed VivaChek about 7 vigilance reports about false positive with different batches of VivaDiag SARS- CoV-2 Ag Rapid Test (REF: VCD05-01-011), most of them regarding batches SE2010037 and SE2011044.

During VivaChek's investigation of all batch records, it was found that the evidence of the implementation of cutting machine cleaning activity which is required to be performed by SOP # after every 3 hours of continuous cutting were missing in the batch records in the manufacturing process batches #SE2010037, #SE2011037, #SE2011044.

Through the test and comparison of 300 clinical samples by retained samples of all batches; one uncompleted T line was identified in each batch of #SE2010037, #SE2011037, #SE2011044. The findings correlate to the complaints received from France and therefore it indicates the violation of cleaning SOP leading to the accumulation of residue on the blade and cutting machine surface. It can further contaminate the T line by the components from C line, results the uncompleted T line which is easy to be read as false "positive".

From our sale record, VivaChek has distributed VivaDiag SARS-CoV-2 Ag Rapid Test (REF:VCD05-01-011) with batch #SE2010037 of 16,050 tests, and batch #SE2011037 of 30,250 tests in Italian market. VivaChek has not distributed batch #SE2011044 in Italian market.

### **Advice on the action to be taken by the user:**

1. For the products, VivaDiag SARS- CoV-2 Ag Rapid Test (REF: VCD05-01-011), of batch #SE2010037 and #SE2011037 sold in Italy, please quarantine immediately, and contact your direct supplier for replacement.
2. Rapid Test is not designed and used for the diagnosis of COVID-19, in most countries, PCR based diagnostic method is the widely accepted as the "Gold Standard" for the confirmation of COVID-19 infection. So, if any doubt on the testing results from Rapid Test, follow-up testing with a molecular diagnostic (PCR) should be considered.